

Vedr.:  Opdatering af software i Mini Crosser X2 brugerpanel	PUP nr.:	25-001
	Udarbejdet af:	LPA
	Dato:	10-02-2025
	Intern <input checked="" type="checkbox"/>	Ekstern <input checked="" type="checkbox"/>

Fra og med serienummer 03-58008-25-02 indføres software revision 33n i X2 brugerpanelet.

Den nye software indlæser default værdierne på max speed, men sikrer også at funktionen er slået fra.

Yderligere muligør den nye software at værdierne kan ændres via kunde specifik USB stikker.

Indtil videre er det kun Norge der har denne funktion slået til. Norges USB fil vil derfor også blive opdateret til ny revision.

<b>MaxSpeedEnabled</b>	<b>= 0; 0=no, 1=yes</b>
MaxSpeedForward	= 15; 6..30 kmh
MaxSpeedReverse	= 5; 3..6 kmh
MaxSpeedKp	= 160; 0..900 units
MaxSpeedKi	= 400; 0..900 units

Software X2 brugerpanel ver 33N kan findes på vores hjemmeside under Tekniker login – Programmering – Programmering Mini Crosser.



**Klaus Smedegaard**  
Administrerende direktør



Regarding:  Updating software in the Mini Crosser X2 user panel	PUP No.:	
	Prepared by:	LPA
	Date:	10-02-2025
	Internal <input checked="" type="checkbox"/>	External <input checked="" type="checkbox"/>

Starting with serial number 03-58008-25-02, software revision 33n is introduced in the X2 user panel.

The new software loads the default values at max speed, but also ensures that the function is switched off.

Furthermore, the new software enables the values to be changed via customer specific USB stick.

So far, only Norway has this function switched on. Norway's USB file will therefore also be updated to the new revision.

<b>MaxSpeedEnabled</b>	= 0; 0=no, 1=yes
MaxSpeedForward	= 15; 6..30 kmh
MaxSpeedReverse	= 5; 3..6 kmh
MaxSpeedKp	= 160; 0..900 units
MaxSpeedKi	= 400; 0..900 units

Software X2 user panel ver 33N can be found on our website under Technician login – Programming – Programming Mini Crosser.



**Klaus Smedegaard**  
General Manager

